


system requirements

Computer requirements

	Planmeca Romexis data base server	Planmeca Romexis client work station
Processor	Pentium 4, 2 Ghz	Pentium 4, 2 GHz
RAM	2 GB	2 GB
Hard disk space	250 GB (RAID)	2 x 80 GB
Graphics card	32 bit graphics card	32 bit graphics card
Monitor	SVGA, 1024 x 900 or more, 24 bit colours or more	SVGA, 1024 x 900 or more, 24 bit colours or more
Network card	Ethernet 10/100 Mbps	Ethernet 10/100 Mbps
Peripherals	CD R/W and/or DVD R/W drive	CD-ROM drive
Backup medium	DAT or other back up device	
Operating system	Windows XP, 2003	Windows XP/2003, Mac OS X*, Linux*
Other	Java platform (Java Virtual Machine 5.0 or later and Java Advanced Imaging)	Java platform (Java Virtual Machine 5.0 or later)

* Mac OS X / Linux support subject to contract

Example installation



*3D imaging
in dental office* 7